

StepMate Repeat Pipettor



Optimized for smooth dispensing of even the smallest volumes

Five easy-to-select settings control dispensing volume, right or left-hand operation

By pushing the slide upwards, the syringe is filled up

The stop lever secures the syringe



Main Data Description

Dial setting	1	2	3	4	5		
Number of doses	48	23	15	11	8		
Syringe volume (ml)	Sample volume (µl)					A %	P %
0.50	10	20	30	40	50	± 0.8	≤ 0.7
1.25	25	50	75	100	125	± 0.8	≤ 0.5
2.50	50	100	150	200	250	± 0.8	≤ 0.4
5.00	100	200	300	400	500	± 0.4	≤ 0.3
12.50	250	500	750	1000	1250	± 0.3	≤ 0.3
25.00	500	1000	1500	2000	2500	± 0.2	≤ 0.3
50.00	1000	2000	3000	4000	5000	± 0.2	≤ 0.2

Ordering information

Cat.No.	Description
72010000	StepMate Repeat Pipettor/Dispenser

EZ-Syringe Tips - Non Sterile

Cat.No.	Description
702371	EZ-Syringe Tips, 0.5ml, 100pk
702373	EZ-Syringe Tips, 1.25ml, 100pk
702375	EZ-Syringe Tips, 2.5ml, 100pk
702377	EZ-Syringe Tips, 5.0ml, 100pk
702379	EZ-Syringe Tips, 12.5ml, 100pk
702381	EZ-Syringe Tips, 25ml, 50pk
702383	EZ-Syringe Tips, 50ml, 25pk

EZ-Syringe Tips - Sterile

Cat.No.	Description
702385	EZ-Syringe Tips, 0.5ml, 100pk
702387	EZ-Syringe Tips, 1.25ml, 100pk
702389	EZ-Syringe Tips, 2.5ml, 100pk
702391	EZ-Syringe Tips, 5.0ml, 100pk
702393	EZ-Syringe Tips, 12.5ml, 100pk
702395	EZ-Syringe Tips, 25ml, 25pk
702397	EZ-Syringe Tips, 50ml, 25pk

E-Z Syringe tips fit all repeating pipettes using typical standard sized syringe tips, such as the Eppendorf® Repeater™ (they do not fit RepeaterPlus™). EZ-Tips are manufactured from virgin polypropylene and high-density polyethylene without lubricants or chemical additives to prevent reagent contamination. Available in non-sterile and sterile. The sterile syringe tips are certified free of DNA, RNase, endotoxins, and ATP to meet the most stringent sample integrity requirements.

Features:

- Lightweight and ergonomic
- One-hand operation
- Factory calibrated
- Volumes from 10µl to 5ml can be dispensed
- Maintenance-free
- Works with seven sizes of disposable, polypropylene syringes from 0.5ml to 50ml
- 3 YEAR WARRANTY

