

CERTIFICATE OF ANALYSIS
No. 36259a

To : S. Schwartz.

Receiving date: 26.05.2020

4 Yehuda H'anachtom St Beer Sheva

Performing date: 04.06.2020

Sample Description	Batch	Manufacturer
SmartSept Hand hygiene solution	B2005133	S. Schwartz


The test was performed according to the Israeli standard -1944 Section: Antimicrobial activity for hand disinfection			
	Performed quality test	Untreated Surface	Treated Surface – Contact time – 1min
			Dilution: NA
1	Pseudomonas aeruginosa 9027	~10 ⁴	0,0,0
2	Staphylococcus aureus –coag.+ 6538	~10 ⁴	1,1,1
3	Escherichia coli 8739	~10 ⁴	1,2,7
4	Streptococcus Faecalis	~10 ⁴	2,0,0
5	Enterobacter Aerogenes 13048	~10 ⁴	1,0,1
6	Salmonella typhimurium 14028	~10 ⁴	3,4,2
7	Aspergillus Spores 16404	~10 ⁴	0,0,0
8	Saccharomyces cerevisiae 2601	~10 ⁴	0,0,0
9	Lactobacillus Plantarum 14917	~10 ⁴	0,0,0
10	Listeria monocytogenes	~10 ⁴	3,2,1

Conclusions: The product is effective against 99.9% of the tested bacteria

- The sampling was done by the customer.
- Terms of maintaining and transporting the sample/s: ambient temperature.
- The laboratory is not authorized to give an opinion and / or interpretation of test results.

Approved by: Dr. S. Sigawi - Laboratory Manager

Date of Approval: 08.06.2020

Signature: 

*****End of the report/ test report *****

- System Laboratories is authorized by the Israel Laboratory Accreditation Authority.
- This certificate is for the costumer use only. Do not use this file, or company name or one of her employees name for commercial or publicity uses without written permission of the lab.
- The report refers only to tested or calibrated subject.
- This report can be copied only in its entire form and only with signed acceptance by the laboratory manager.