

# CERTIFICATE OF QUALITY CONTROL

<b>Description</b>	: Water Sampling Bottles - with sodium thiosulfate -
<b>Catalog Number</b>	: C.041.80.504
<b>Lot Number</b>	: BWPPS21901301
<b>Manufacture Date</b>	: April 2021
<b>Expiration Date</b>	: April 2024
<b>Material</b>	: Polypropylene
<b>Color</b>	: Clear
<b>Leakage Test</b>	: Pass
<b>Dimensional Control</b>	: Pass
<b>Leakage Test</b>	: Pass
<b>Visual Appearance</b>	: Within acceptability limits
<b>Sterilization</b>	: Gamma irradiation according to 11137-1:2015 and 13485:2016 standards
<b>Dimensional Control</b>	: Pass
<b>Packing Type</b>	: Bulk packed

Catalog Number	Volume	w X d X H	Neck Inner Diameter	Neck Size
	ml	mm	mm	
C.041.80.504	1000	90 X 90 X 175	39	GL 50

We declare that the above mentioned batch on this certificate is manufactured under the ISO 9001:2015 Quality Management System and sent to gamma irradiation for sterilization. Each bottle includes 20 mg/lit sodium thiosulfate according to ISO/FDIS 19458:2006. Bottles are individually labelled specifying e.g. sample and source descriptions, code and lot numbers and expiry date.

The above mentioned properties tested ISOLAB Quality Control Laboratory.  
 For and on behalf of

**Quality Department**  
**ISOLAB LABORGERÄTE GmbH**  
 Date: April 2021

**Quality Manager**

A handwritten signature in black ink, appearing to be 'R. Schmitt', located in the upper right corner of the page.

Am Dillhof 2 - 63863 Eschau / GERMANY

Tel : +49 93 74 / 978 55-0 Fax : +49 93 74 / 978 55-29 Web: [www.isolab.de](http://www.isolab.de) E-mail: [info@isolab.de](mailto:info@isolab.de)