

FICHE TECHNIQUE

GRADE PE45

DESIGNATION

- Papier de pesée
- Matière : 100% fibres de cellulose
- Résistant à l'humidité
- Liant : aucun

APPLICATIONS

- Pesée précise de toutes sortes de matières (liquides, poudres...)

SPECIFICATIONS TECHNIQUES

CARACTERISTIQUES	VALEURS	UNITES	NORMES
Grammage	45	g/m ²	DIN 53104
Epaisseur	0.06	mm	
Résistance à l'éclatement	0.8	kg/cm ²	

DONNEES COMMERCIALES

- Disponibles en feuilles

Dimensions (mm)	Conditionnement	Réf. commerciale
50x50	250	PE45F0505
76x76	500	PE45F0707D
95x110	250	PE45F0911
95x110	100	PE45F0911C
100x100	250	PE45F1010
102x102	500	PE45F1010D
100x100	1000	PE45F1010M
120x120	250	PE45F1212
150x150	250	PE45F1515
152x152	500	PE45F1515D

TECHNICAL DATA

GRADE PE45

DESCRIPTION

- Weighing paper
- Material: 100% cellulose
- Moisture resistant
- Binder: no

APPLICATIONS

- Precise weighing of all kinds of materials (liquids, powders, etc.)

TECHNICAL SPECIFICATIONS

CHARACTERISTICS	VALUES	UNITS	STANDARDS
Weight	45	g/m ²	DIN 53104
Thickness	0.06	mm	
Burst strength	0.8	kg/cm ²	

SALES DATA

- Available in sheets (F)

Size (mm)	Packing	Commercial ref.
50x50	250	PE45F0505
76x76	500	PE45F0707D
95x110	250	PE45F0911
95x110	100	PE45F0911C
100x100	250	PE45F1010
102x102	500	PE45F1010D
100x100	1000	PE45F1010M
120x120	250	PE45F1212
150x150	250	PE45F1515
152x152	500	PE45F1515D