

FICHE TECHNIQUE

REF PH114PS

DESIGNATION

- Bandelettes plastiques
- pH 0.0-14.0

CARACTERISTIQUES

- Support : bandelette plastique + papier imprégné à base de réactifs
- Aspect : solide (bandelette plastique)
- Non contaminant (pas de migration de l'échantillon par des composants)

CONSERVATION

- Au sec et à l'abri de la lumière

APPLICATIONS

- Analyses de solutions en milieu faiblement acide, neutre ou légèrement alcalin

PRECISION DE LA MESURE

- La mesure est étalonnée par saut de 1.0 unité de pH
- Le papier réactif réagit par méthode colorimétrique sur les 4 zones indicatrices en comparant le résultat obtenu à l'échelle fournie.

DONNEES COMMERCIALES

- Disponible en boîtes de 100 ou 50 bandelettes plastiques

TECHNICAL DATA

REF PH114PS

DESCRIPTION

- Indicator plastic strips
- pH0.0-14.0

CHARACTERISTICS

- Material: plastic strip + reagent impregnated filter paper
- Appearance: Solid (plastic strip)
- Non-bleeding (no contamination of the tested sample)

STORAGE

- Dry and light-free storage

USE

- Analysis of solutions in weakly acidic, neutral or slightly alkaline medium

ACCURACY OF MEASUREMENT

- The measurement is calibrated every 1.0 pH unit.
- The reagent paper reacts by colorimetric method on 4 coloured zones; the results obtained must be compared to the colorimetric scale given in each box

SALES INFORMATION

- Available in boxes of 100 or 50 plastic strips