

## YP Series Electronic Balance



### Product introduction:

YP series electronic balance has fluent and novel appearance design which reflects the pursuit for modern concept. New design brings faster weighing response, higher stability and stronger anti-interference performance. Our YP series electronic balance can be customized upon user's request. With strong flexibility, it will be the best choice for field operation.

### Feature :

- AC/DC exchangeable, balance level, backlight LCD display, rapid weighing response and high stability.
- Sensor adopts unique locking device, which guarantees protection of the balance during transportation.

### Function:

- System operations: including auto-counting, percentage, automatic calibration, multi-point calibration, standard weights (optional).
- Multiple functions guarantee easy and reliable operation, including unit conversion (metric, carat, gold ounce, etc.), stability, memory, full-scale ranging.
- Built-in RS232C output interface guarantees direct connection to external equipment such as computer and printer.

Model	Weighing Range	Readable Precision	Outline Dimension	Outline Dimension	Scale Size	Net Weight
YP1002	0-100g	10mg	250*180*60mm	365*265*170mm	Ø125mm	1kg
YP2002	0-200g	10mg	250*180*60mm	365*265*170mm	Ø125mm	1kg
YP3002	0-300g	10mg	250*180*60mm	365*265*170mm	Ø125mm	1kg
YP4002	0-400g	10mg	250*180*60mm	365*265*170mm	Ø125mm	1kg
YP5002	0-500g	10mg	250*180*60mm	365*265*170mm	Ø125mm	1kg
YP6002	0-600g	10mg	250*180*60mm	365*265*170mm	Ø125mm	1kg
YP10002	0-1000g	10mg	250*180*60mm	365*265*170mm	130*130mm	1kg
YP20002	0-2000g	10mg	250*180*60mm	365*265*170mm	130*130mm	1kg
YP30002	0-3000g	10mg	250*180*60mm	365*265*170mm	160*160mm	1kg
YP50002	0-5000g	10mg	250*180*60mm	365*265*170mm	160*160mm	1kg
YP2001	0-200g	100mg	250*180*60mm	365*265*170mm	Ø125mm	1kg
YP4001	0-400g	100mg	250*180*60mm	365*265*170mm	Ø125mm	1kg
YP5001	0-500g	100mg	250*180*60mm	365*265*170mm	Ø125mm	1kg
YP6001	0-600g	100mg	250*180*60mm	365*265*170mm	Ø125mm	1kg
YP8001	0-800g	100mg	250*180*60mm	365*265*170mm	Ø125mm	1kg
YP10001	0-1000g	100mg	250*180*60mm	365*265*170mm	130*130mm	1kg
YP16001	0-1600g	100mg	250*180*60mm	365*265*170mm	130*130mm	1kg
YP20001	0-2000g	100mg	250*180*60mm	365*265*170mm	130*130mm	1kg
YP30001	0-3000g	100mg	250*180*60mm	365*265*170mm	160*160mm	1.3kg
YP40001	0-4000g	100mg	250*180*60mm	365*265*170mm	160*160mm	1.3kg
YP50001	0-5000g	100mg	250*180*60mm	365*265*170mm	160*160mm	1.3kg
YP60001	0-6000g	100mg	250*180*60mm	365*265*170mm	160*160mm	1.3kg
YP80001	0-8000g	100mg	250*180*60mm	365*265*170mm	160*160mm	1.3kg
YP1000	0-1000g	1g	250*180*60mm	365*265*170mm	130*130mm	1kg
YP2000	0-2000g	1g	250*180*60mm	365*265*170mm	130*130mm	1kg
YP3000	0-3000g	1g	250*180*60mm	365*265*170mm	160*160mm	1kg
YP6000	0-6000g	1g	300*120*160mm	365*265*170mm	160*160mm	1kg