

Desktop Drying Oven And Incubator

Features

1. 18L small desktop drying oven/incubator.
2. Stainless steel corrosion-resistant inner chamber.
3. Intelligent Temperature control instrument adopts PID control program, dual digital display,
4. With timer (0-9999min) and over-temperature alarm.
5. Heat isolation door with tempered glass observation window.



Specifications

Model		WHL-25AB	WP-25AB
Heating type		neutral convection	
Performance	temp. range	RT+10~300 °C	RT+5~65 °C
	temp. resolution	0.1 °C	
	temp. fluctuation	±0.1°C	
	temp. accuracy	±3.5%	
Structure	chamber material	stainless steel	
	external material	cold rolled steel, powder coating	
	thermal insulation	aluminium silicate fiber	
	heater	stainless steel heating tube	
	rated power	0.6KW	0.18KW
	air outlet	inner diameter 28mm*1, top	
Controller	temp. control	fuzzy PID control	
	temp. display	LED display	
	timer	0-9999 minutes	
	sensor	Pt100	
Safety device		over-temperature alarm	
Specification	inner (W*D*H) / mm	260*260*260	
	exterior (W*D*H) / mm	410*410*620	
	packing (W*D*H) / mm	530*560*700	
	capacity	18L	
	shelf load-bearing	15KG	
	shelf number	2 (standard) 6 (max)	
	voltage	AC220V / 2.7A	AC220V / 0.8A
	N.W / G.W (kg)	27/35	
Optional configuration		shelf	