

## GW2145

---



LABORATORY GLASSWARE WASHER  
BASIC LINE - GLASSWARE WASHER 45 CM  
[Download PDF](#)

---

- ✓ washing up to 95°C - antibacterial thermal disinfection
- ✓ programs: 9 default + 6 custom
- ✓ microprocessor programmable electronic control
- ✓ serial port RS232 for connection to external printer or PC
- ✓ active thermodynamic drying
- ✓ alkaline powder detergent dispenser
- ✓ peristaltic pump for acid liquid neutralizer
- ✓ possibility of installing 1 additional peristaltic pump
- ✓ washing chamber and internal door in AISI 316L stainless steel
- ✓ external panels made of AISI 304 stainless steel
- ✓ n. 2 independent washing levels through telescopic guides
- ✓ external dimensions LxDxH: 450 x 620 x 850 mm
- ✓ internal working dimensions LxDxH: 370 x 480 x 570 mm
- ✓ incorporated water softener with salt sensor
- ✓ cold water and demineralized water connections
- ✓ electronic door lock system for user safety
- ✓ electrical connection: 1/N/PE 230 V - 50 Hz - 3,3 kW max
- ✓ 60 Hz version available
- ✓ Schuko plug 16 A - 250 V included
- ✓ upper and lower trolleys D-CS1 and D-CS2 included

Smeg S.p.A. Via Leonardo Da Vinci, 4 - 42016 Guastalla (RE) Italy - P. IVA 01555030350  
Phone +39 0522 8211 - fax +39 0522 821453 - [service.instruments@smeg.it](mailto:service.instruments@smeg.it)

